

LOCKING PHILOS PLATE


3.5MM

- Anatomically contoured design for accurate fit
- Multiple locking screw options for enhanced angular stability
- Ideal for osteoporotic and complex fractures
- Low-profile plate to minimize soft-tissue irritation



ISO 9001:2015
ISO 13485:2016
WHO GMP
CE Certified

TRAUMA | SPINE | CMF




Locking Philos Plate 3.5mm

Cat No: (S.S.) S00.010.003 to 014 | (T.T.) T00.010.003 to 014

Size / Holes	03, 04, 05, 06, 07, 08, 09, 10, 11, 12, 13, 14 Holes
Material	Stainless Steel (SS 316L) / Titanium
Brand	OnnixMed
Indication	This plate is used for stabilizing proximal humerus fractures, especially complex, dislocated, or osteoporotic ones, providing stable, fixed-angle support for early mobilization. Its design allows for angled screws, securing bone fragments even in poor bone
Used With	Cortical Screw 3.5mm, Locking Screw 3.5mm, Locking Cancellous Screw 3.5mm, Locking Cancellous Screw 3.5mm and 4.0mm

Cortical Screw Ø 3.5mm

Cat No: (S.S.) S02.008.010 to 100 | (T.T.) T02.008.010 to 100



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

Locking Cortical Screw Ø 3.5mm

Cat No: (S.S.) S02.028.010 to 100 | (T.T.) T02.028.010 to 100




Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

Locking Cancellous Screw Ø 3.5 / Ø 4.0mm

Cat No: (S.S.) S02.029.010 to 090 | (T.T.) T02.029.010 to 090

Cat No: (S.S.) S02.030.010 to 090 | (T.T.) T02.030.010 to 090



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 35, 36, 38, 40, 42, 44, 45, 46, 48, 50, 52, 54, 55, 56, 58, 60, 65, 70, 75, 80, 85, 90 mm

CONTACT US

ONNIX TechnoMedicals

- Registered Office: No.5, Street No. 15A, Block A2, West Sant Nagar, Burari, Delhi - 110084
- Corporate Office: Kh. No. 22/27/1 & 22/27/2, GF & FF, Street No. 09, Near Heater Factory, Samaypur Industrial Area, Delhi - 110042

- Call us: +91-7701851136, +91-8800263670, +91-8851719001
- Email: info@onnixmed.com, onnixtechmed@gmail.com, sales@onnixmed.com
- Website: www.onnixmed.com