

# LOCKING CLAVICLE HOOK PLATE - 3.5/4.0MM (LEFT & RIGHT)

- Anatomically contoured design
- Superior subacromial support
- Locking screw technology for enhanced stability
- Helps maintain reduction during healing
- Available in Medical grade Stainless Steel 316L / Titanium



ISO 9001:2015  
ISO 13485:2016  
WHO GMP  
CE Certified

**TRAUMA | SPINE | CMF**



Dedicated to Quality. Committed to Service.



## Locking Clavicle Hook Plate - 3.5/4.0mm (LEFT & RIGHT)


**Cat No (LEFT):** (S.S.) S00.006.003 to 009 | (T.T.) T00.006.003 to 009

**Cat No (RIGHT):** (S.S.) S00.007.003 to 009 | (T.T.) T00.007.003 to 009

Size / Holes	03, 04, 05, 06, 07, 08, 09 Holes
Material	Stainless Steel (SS 316L) / Titanium
Brand	OnnixMed
Indication	This plate is used to fix severe clavicle fractures, especially at the outer (lateral) end, often with acromioclavicular (AC) joint disruption, providing stable fixation by combining locking screws for angular stability with a hook to help restore the joint's height and stability, minimizing soft tissue irritation.
Used With	Cortical Screw 3.5mm, Locking Screw 3.5mm

### Cortical Screw Ø 3.5mm


**Cat No:** (S.S.) S02.008.010 to 100 | (T.T.) T02.008.010 to 100



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

### Locking Cortical Screw Ø 3.5mm

**Cat No:** (S.S.) S02.028.010 to 100 | (T.T.) T02.028.010 to 100



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

### Locking Cancellous Screw Ø 3.5mm

**Cat No:** (S.S.) S02.029.010 to 090 | (T.T.) T02.029.010 to 090



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 35, 36, 38, 40, 42, 44, 45, 46, 48, 50, 52, 54, 55, 56, 58, 60, 65, 70, 75, 80, 85, 90 mm

## CONTACT US

### ONNIX TechnoMedicals

- Registered Office: No.5, Street No. 15A, Block A2, West Sant Nagar, Burari, Delhi - 110084
- Corporate Office: Kh. No. 22/27/1 & 22/27/2, GF & FF, Street No. 09, Near Heater Factory, Samaypur Industrial Area, Delhi - 110042

- Call us: +91-7701851136, +91-8800263670, +91-8851719001
- Email: [info@onnixmed.com](mailto:info@onnixmed.com), [onnixtechmed@gmail.com](mailto:onnixtechmed@gmail.com), [sales@onnixmed.com](mailto:sales@onnixmed.com)
- Website: [www.onnixmed.com](http://www.onnixmed.com)