

LOCKING POSTEROMEDIAL PROXIMAL TIBIAL PLATE 3.5MM

- Anatomical posteromedial contouring for optimal fit
- Locking screw technology for enhanced stability
- Supports complex & high-stress fracture patterns
- Minimizes soft tissue irritation
- Trusted performance for consistent surgical outcomes
- Available in medical grade SS 316L / Titanium



ISO 9001:2015
ISO 13485:2016
WHO GMP
CE Certified

TRAUMA | SPINE | CMF



Dedicated to Quality. Committed to Service.




Locking Posteromedial Proximal Tibial Plate 3.5mm

Cat No: (S.S.) S00.065.001 to 010 | (T.T.) T00.065.001 to 010

Size / Holes	01, 02, 03, 04, 05, 06, 07, 08, 09, 10 Holes
Material	Stainless Steel (SS 316L) / Titanium
Brand	OnnixMed
Indication	This plate is used for fixation of posteromedial proximal tibia fractures including buttressing of fractures of the proximal, distal, and metaphyseal areas of the tibia.
Used With	Cortical Screw 3.5mm, Locking Screw 3.5mm, Locking Cancellous Screw 3.5mm and 4.0mm

Cortical Screw Ø 3.5mm


Cat No: (S.S.) S02.008.010 to 100 | (T.T.) T02.008.010 to 100



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

Locking Cortical Screw Ø 3.5mm

Cat No: (S.S.) S02.028.010 to 100 | (T.T.) T02.028.010 to 100



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65, 70, 75, 80, 85, 90, 95, 100 mm

Locking Cancellous Screw Ø 3.5 / Ø 4.0mm

Cat No: (S.S.) S02.029.010 to 090 | (T.T.) T02.029.010 to 090

Cat No: (S.S.) S02.030.010 to 090 | (T.T.) T02.030.010 to 090



Material	Stainless Steel (SS 316L) / Titanium
Size (Length)	10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 35, 36, 38, 40, 42, 44, 45, 46, 48, 50, 52, 54, 55, 56, 58, 60, 65, 70, 75, 80, 85, 90 mm

CONTACT US

ONNIX TechnoMedicals

- Registered Office: No.5, Street No. 15A, Block A2, West Sant Nagar, Burari, Delhi - 110084
- Corporate Office: Kh. No. 22/27/1 & 22/27/2, GF & FF, Street No. 09, Near Heater Factory, Samaypur Industrial Area, Delhi - 110042

- Call us: +91-7701851136, +91-8800263670, +91-8851719001
- Email: info@onnixmed.com, onnixtechmed@gmail.com, sales@onnixmed.com
- Website: www.onnixmed.com